**We have solutions for your research purposes**

Water Maze

CPP

Open Field

Plus Maze

T-Maze

Rotarod

Analgesia Hotplate


Grip Strength Meter

Rotameter

Electronic von Frey

Automatic Shuttle Box

Wheel Trainer

Programmable  
Automatic Treadmill  
1,5,10 Channel

**TECHNIC AZMA** is Producer Of Test Equipment In Field Of Behavioral And Neuroscience For Physiology And Pharmacology Studies With **MazeRouter**'s Name

Address : Floor 3 / No.40 / Farabi Ave. / Tabriz / IRAN

Telfax : 041- 36556830 Mobile : 0914 419 2012

## Automatic Small Animal Treadmill

### General

MazeRouter provides a complete solution enabling forced exercise training and accurate testing of fatigue in rodents. Different treadmill models are available depending on the user's needs (manual and Automatic in 1/5/10 lane).

Basically, the treadmill apparatus consists of a rolling belt with adjustable speed (2 ~ 80 m/min) and slope ( from -20 to 20 degrees) and a control unit with touch screen and large graphical lcd display.

The treadmill unit controls the speed of the belt, shows measured data in its display ( belt speed, distance travelled , number of shocks and accumulated duration) and provides adjustable electrical shock to the grid.

The electrical shock supplied by the grid is of constant intensity (from 0.1 to 9.9 mA), that is, the current which circulates through the animal (and therefore its effect) only depends on the value of the mA chosen and not of the subject (quantity of body mass in contact with the bars, perspiration, etc.).

The rolling belt is built with especially selected materials to guarantee the best performance under conditions of intensive use and the minimum operations of maintenance.

The lanes have sufficient width for the subject to correct its errors in coordination, thereby allowing an exact measurement of the fatigue without deficiencies in motor coordination.

### Main Features:

- Touch screen micro-processor controller
- Plexiglass restrainer
- Rolling belts with adjustable speed and slope
- Silent operation, even at high regimes
- Electrical adjustable slope of rolling belt
- Accurate control of the intensity of the shock delivered
- Positive and negative slope
- Easy to clean

### Applications:

Improvement of sportive human performances, Exercise, Fatigue, Oxidative stress, Respiratory metabolism (O<sub>2</sub>/CO<sub>2</sub>), Diabetes, Parkinson, Ischemia, Osteopenia/osteoporosis...

## Shuttle box (Active/Passive Avoidance)

### General

shuttle box system is fear-motivated tests classically used to assess short-term or long-term memory on small laboratory animals (rat, mice).


Basically, active/passive avoidance working protocols involve timing of transitions, i.e. time that the animal takes to move from the white compartment to the black one after a conditioning session -in which the entry into the black compartment is punished with a mild inescapable electrical shock- is carried out.

MazeRouter Shuttle Box controlled through SB100 Micro-Processor Base Controller with touch screen.

### Main Features

- Latency to enter into the black compartment is the parameter that measure
- Touch screen micro-processor controller
- silent and Automatic guillotine gate
- Specific models for rats or mice
- Removable stainless steel shock grid
- best design for simple cleaning
- A shocker with adjustable current generator

The MazeRouter Shuttle Box System is defined by a white illuminated compartment and a black dark compartment separated by a guillotine gate. The animal's position is detected by using high sensitivity photo Electric transducers providing higher effective and reliable detection of animal responses (zones entries).

**Removable cage from shock grid and excrement tray**

**IDEAL TO STUDY  
Short-term or  
Long-term memory**

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## Hot Plate Analgesia Meter

### General

The MazeRouter hotplate Analgesia Meter has been designed to perform rapid precise screening of narcotic type analgesic drugs (morphine, codeine, etc.) on the animal by the hot plate test.

This method, devised by N.B. Eddy and D. Leimbach evaluates the reaction time of mice which are dropped on a heated surface and thus confronted with a heat stimulus applied to the plantar surface.

When a central analgesic is administered to the animals this reaction time is markedly increased.

### Main Features

- ▶ Touch screen micro-processor controller
- ▶ Plexiglass
- ▶ Wireless Remote control to control the experiment from far away
- ▶ Removable hot platform for easy cleaning the surface of plate
- ▶ Temperature control from 20 ~ 80°C with Holding accuracy +/- 1°C
- ▶ Specific models for mice and rat

### Hot Plate Analgesia Meter

The time taken between contacting the heated measuring platform and the appearance of the typical "paw licking reaction" is a measure of the sensitivity of the animal to pain .(plantar test)

The application of the analgesically effective substances increases the latency time by increasing the pain threshold compared with untreated animals.

## Small Animal Rotarod (mice and Rat )

### General

The MazeRoter Rotarod is used to assess the effect of drugs on the motor coordination or fatigue resistance on mice and rats using one unit in a safe and humane way. The animals are placed on textured drums to avoid slipping . When an animal drops onto the individual sensing platforms below, test results are recorded.

when a rat (or mic) falls off its cylinder section on to the plate be-low, the plate trips thereby recording the animal's endu-rance time in seconds . A large, very readable backlit graphic display shows the actual angular speed (RPM). At the end of a run, the display shows for each animal the running time and the instrument rotation speed at the time that animal fell off.

The RotaRod is controlled by an advanced microprocessor which provides precise timing control and accurate speed regulation , speed range is 2 ~ 45 RPM. Rotation can be electronically set at a constant speed by toch screen display on the front panel . Alternatively , acceleration rate may be selected at a defined time .

### Main Features

- Very stable with metal framework
- Touch screen micro-processor controller
- Silent operation, even at high regimes
- 1 ~ 45 RPM
- accurate speed regulation
- selectable acceleration rate
- Micro-switch detection of fall
- Individual lane timers
- Automatic recording of latencies to fall and rotation speed

### Applications:

Aging  
Fatigue  
Exercise  
Motor activity  
Oxidative stress  
Drug screening  
Parkinson disease  
Huntington disease  
Alcohol dependence  
Improvement of sportive performances  
Respiratory metabolism (O<sub>2</sub>/CO<sub>2</sub>)  
Diabetes ,  
Ischemia ,  
Osteopenia / osteoporosis...

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Frequency = 50 Hz  
Duty cycle = 50% , 75% , 25%  
Voltage = 5 v  
Current = 5mA  
Number of pulse = 3  
Delay time = 1.5 sec  
Number of trial = 3

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- Number of Lick Before Shock (1~50)
- Shock Time (0.1 ~ 5.0 Sec)
- Shock Current (0.1 ~ 5.0 mA)
- Recording Interval (1~60 min)
- Experiment Duration (60 ~ 900 min)
- Start Mode (Start With Run Key / Start With First Shock / Start With First Lick)

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## Grip Strength Test

helps determine the maximal peak force developed by a rodent (rat or mouse)

### General

To assess muscle hyperalgesia this instrument measures the forelimb grip force using a digital force transducer.

The Grip Strength test is a perfect complement to the gold standard MazeRouter Rota-Rod device for motor coordination and motor function experiments.

The effects of drugs, toxins, muscle relaxants, disease, ageing or neural damage on muscle strength may be assessed.

The transducer is connected to a wire mesh grid connected to a base plate. The animal is held by its tail and is gently passed over the mesh until it grasps the grid with its front paws.

Three grip force measures are made and the average of these readings represents the animal's forelimb grip force at that particular time. The last maximum force is held in a "peak & hold" type readout until reset.

The unit gives readings in one gram increments. Both rat and mouse can make Test with this instrument. Maximum force range is 2000 Gr (50% over range allowed). **Higher ranges available on special order.**

The heavy base plate has suction feet to resist even large pulling forces.

### Main Features

- ▶ Easy to use
- ▶ No calibration needed
- ▶ Automatic Peak Detector
- ▶ 0 ~ 2000 gr in 1gr resolution
- ▶ Stand alone system, no PC needed
- ▶ Touch screen micro-processor controller

### Applications:

Aging  
Neuromuscular diseases  
Phenotyping  
Drug screening  
Parkinson disease  
Huntington disease

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
- Subject name
- Current zone
- Current speed
- Maximum speed
- Distance travelled
- Time in zone
- Number of visit
- Group name
- Group number
- Day
- Animal number
- Trial
- Time remaining

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

## Conditioned Place Preference Box (CPP)

### General

The new MazeRouter Conditioned Place Preference (CPP) is a 3-compartment box to evaluate the abuse potential of substances and the motivational effects of drugs. The box includes the contextual cues required by the experimental paradigm.

The 3 compartments differ for the walls colour and patterns and for the floor patterns and texture. The floors are interchangeable so that the tactile difference between the 3 compartments can be easily adjusted by the scientist.

The new CPP box has been designed and optimized for use with MazeRouter video-tracking software or manual scoring.

### Main Features

- ▶ Interchangeable patterned floors
- ▶ Striped and dark compartment
- ▶ Optimized for Video-tracking
- ▶ Specific models for rats or mice
- ▶ Designed for multiple-cage systems

### Conditioning Place Preference (CPP)

Center Arm Box

Side Arm Box

**Interchangeable  
floors for tactile  
stimulation**

**IDEAL TO STUDY  
Drug Abuse  
Addiction**

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